

# Autoclave Training



**URSINUS COLLEGE  
BIOLOGY AND CHEMISTRY  
DEPARTMENTS**

# Objectives



- What materials can be autoclaved?
- What Personal Protective Equipment (PPE) is necessary for autoclave use?
- How do you properly load and unload the autoclave?
- Who do you go to for assistance?

# Importance of Proper Autoclave Use



- Tremendous pressure from steam in the chamber provides explosive potential.
- High temperatures (250°F) and presence of extremely hot water creates potential for burns and scalding.
- Inadequate decontamination allows for the potential of biological hazards to contaminate personnel and the environment.


# Autoclave Explosion



# What can be autoclaved ?




- Glassware
- Polypropylene plastic tubes and pipette tips
- Aqueous Solutions and water
- Non-Biohazardous laboratory waste



# Principles of Autoclave Operation



- Steam penetrates objects in the autoclave
- Condensation creates negative pressure and draws in additional steam
- Moist heat kills microorganisms via coagulation of proteins



# Autoclave Use Basics




- Personal Protective Equipment (PPE)
- Loading
- Operating
- Unloading
- Improper Autoclave Practices



# PPE for Autoclave Users


- Eye Protection
- Buttoned lab coat is strongly recommended
- Closed-toed Shoes
- Heat-resistant Gloves



# Prior to Loading the Autoclave



- Be sure the material can be autoclaved
  - No corrosives (acids, bases, phenol)
  - No solvents or volatile chemicals (ethanol, methanol, chloroform, bleach)
  - No sharps or red bag waste
- Choose appropriate containers and autoclave bags
  - Clear bags from inside the cardboard receptacles
  - Do not overfill bags or containers
  - Do not allow the plastic to touch the sides of the autoclave



# Loading the Autoclave



- To allow for steam penetration
  - Loosely close autoclave bags
  - Loosen lids on bottles containing solutions to prevent the shattering of glass during pressurization
- Place items in an autoclavable polypropylene bin
  - Bottles should not be filled more than 2/3 & keep 1-2 inches between bottles
  - Cover open tops of glassware with foil
  - Label with heat sensitive indicator tape



# Loading the Autoclave



- Always autoclave clean items and waste separately.
- Note the cycle end time and come back to collect your items promptly. Close the autoclave door and turn the autoclave OFF when you are finished.
- Items left unlabeled for an excessive period of time will probably disappear!



# Operating the Autoclaves



- Be sure the autoclave is functioning properly before use
- Close door properly and securely
- Choose the correct cycle for your material

## **Thomas Hall autoclave cycles**

- P07, gravity 15/15 for glassware
- P08, liquid 15/15
- P09, liquid 30/30
- P10, liquid 45/45
- P11, gravity 30/15 for lab waste

# Operating the autoclaves



## Pfahler Autoclave cycles


Program	Temp	Exhaust	Dry
P01	250F	30 min	30 min Standard dry materials (glassware/plastics)
P02	250F	10 min	5 min
P03	275F	3 min	0 min
P04	275F	10 min	0 min
P05	250F	20 min	Liquid Standard liquid run
P06	250F	15 min	Liquid

# Unloading the Autoclave



- Wear the recommended Personal Protective Equipment
- Allow the autoclave to completely finish cycle
  - Pressure gauge must read zero
- Open door slightly to allow steam to escape
- Carefully remove items
  - Be especially careful with liquids
  - Transport hot liquids in secondary containment
- Be considerate and return autoclave bins and gloves

# Improper Autoclave Practices



Flip flops



Biohazards



Ethanol

# Results of Improper Autoclave Use



# Review



- What materials can be autoclaved ?

**Glassware ( beakers, graduated cylinders, flasks etc.)**

**Polypropylene plastic tubes, and pipette tips**

**Aqueous solutions and water**

**Non-Biohazardous Laboratory Waste (From inside the cardboard waste receptacles, may contain petri dishes, agar, bacteria, *C.eligans*)**

# Review

- What Personal Protective Equipment (PPE) is necessary for safe autoclave use?



## Review



- How do you properly load the autoclave?

**Load material to allow efficient steam penetration.**

**Autoclave clean items and waste separately.**

**Do not allow material to be autoclaved to touch the sides or top of the chamber.**

**Use secondary containment.**

# Review



How do you properly unload the autoclave?

**Wear safety glasses, lab coat, closed-toed shoes, orange autoclave gloves.**

**Retrieve your items on time.**

**Transport the hot liquids in secondary containment.**


# Who to ask for assistance?



- Your faculty mentor
- Ann Breen, Biology Laboratory Manager x3072
- Brian Phillips, Chemistry Laboratory Manager x3346

# Summary

- Contact your faculty mentor to schedule your hands-on autoclave training and quiz.
- You can operate the autoclave following hands-on training, a review of this PowerPoint presentation and passing the quiz.
- Your quiz will be kept on file in the Biology or Chemistry Dept.



# References



University of Kentucky, Dept. of Biological Safety

University of South Carolina Autoclave Safety Policy

UNC Environmental Health and Safety Dept.